The logo for amoena, featuring the brand name in a lowercase, sans-serif font. The background of the top right corner is a teal-to-white gradient triangle pointing towards the top right.

amoena

Redonner Confiance

The title 'LIVRE BLANC MÉDICAL' is displayed in a large, teal, sans-serif font, centered on the page. The background is white with teal geometric shapes in the corners.

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5 cf. Laura, S., Clark, D. and Harvey, F. (2004), Patient preference for bra or binder after breast surgery. ANZ Journal of Surgery, 74: 463-464. <https://doi.org/10.1111/j.1445-2197.2004.03019.x>

6 cf. Nicklaus, K.M., Bravo, K., Liu, C. et al. Undergarment needs after breast cancer surgery: a key survivorship consideration. Support Care Cancer 28, 3481–3484 (2020). <https://doi.org/10.1007/s00520-020-05414-z>

7 cf. Hansdorfer-Korzon R, Teodorczyk J, Gruszecka A, Wydra J, Lass P. Relevance of low-pressure compression corsets in physiotherapeutic treatment of patients after mastectomy and lymphadenectomy. Patient Prefer Adherence. 2016;10:1177-1187 <https://doi.org/10.2147/PPA.S108326>

8 Lasinski B. B. (2013). Complete decongestive therapy for treatment of lymphedema. Seminars in oncology nursing, 29(1), 20–27. <https://doi.org/10.1016/j.soncn.2012.11.004>

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2 cf. Nicklaus, K.M., Bravo, K., Liu, C. et al. Undergarment needs after breast cancer surgery: a key survivorship consideration. Support Care Cancer 28, 3481–3484 (2020). <https://doi.org/10.1007/s00520-020-05414->

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1 cf. Todd M. (2017). Identification, assessment and management of breast oedema after treatment for cancer. International journal of palliative nursing, 23(9), 440-444. <https://doi.org/10.12968/ijpn.2017.23.9.440>

2 cf. Boughey, J. C., Hoskin, T. L., Cheville, A. L., Miller, J., Loprinzi, M. D., Thomsen, K. M., Maloney, S., Baddour, L. M., & Degnim, A. C. (2014). Risk factors associated with breast lymphedema. Annals of surgical oncology, 21(4), 1202-1208. <https://doi.org/10.1245/s10434-013-3408-5>

3 cf. DiSipio, T., Rye, S., Newman, B., & Hayes, S. (2013). Incidence of unilateral arm lymphoedema after breast cancer: a systematic review and meta-analysis. The Lancet. Oncology, 14(6), 500-515. [https://doi.org/10.1016/S1470-2045\(13\)70076-7](https://doi.org/10.1016/S1470-2045(13)70076-7)

4 cf. Abouelazayem, M., Elkorety, M., & Monib, S. (2021). Breast Lymphedema After Conservative Breast Surgery: An Up-to-date Systematic Review. Clinical breast cancer, 21(3), 156-161. <https://doi.org/10.1016/j.clbc.2020.11.017>

5 cf. Abouelazayem, M., Elkorety, M., & Monib, S. (2021). Breast Lymphedema After Conservative Breast Surgery: An Up-to-date Systematic Review. Clinical breast cancer, 21(3), 156-161. <https://doi.org/10.1016/j.clbc.2020.11.017>

6 cf. Verbelen, H., Tjalma, W., Dombrecht, D. et al. Breast edema, from diagnosis to treatment: state of the art. Arch Physiother 11, 8 (2021). <https://doi.org/10.1186/s40945-021-00103-4>

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8 cf. Hansdorfer-Korzon R, Teodorczyk J, Gruszecka A, Wydra J, Lass P. Relevance of low-pressure compression corsets in physiotherapeutic treatment of patients after mastectomy and lymphadenectomy. Patient Prefer Adherence. 2016;10:1177-1187. <https://doi.org/10.2147/PPA.S108326>

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9 cf. Hansdorfer-Korzon R, Teodorczyk J, Gruszecka A, Wydra J, Lass P. Relevance of low-pressure compression corsets in physiotherapeutic treatment of patients after mastectomy and lymphadenectomy. Patient Prefer Adherence. 2016;10:1177-1187 <https://doi.org/10.2147/PPA.S108326>

10 cf. Verbelen, H., Tjalma, W., Dombrecht, D. et al. Breast edema, from diagnosis to treatment: state of the art. Arch Physiother 11, 8 (2021). <https://doi.org/10.1186/s40945-021-00103-4>

11 cf. Executive Committee of the International Society of Lymphology (2020). The diagnosis and treatment of peripheral lymphedema: 2020 Consensus Document of the International Society of Lymphology. Lymphology, 53(1), 3-19.

12 cf. Gregorowitsch, M. L., Van den Bongard, D., Batenburg, M., Traa-van de Grootevheen, M., Fuhler, N., van Het Westeinde, T., van der Pol, C. C., Young-Afat, D. A., & Verkooijen, H. M. (2020). Compression Vest Treatment for Symptomatic Breast Edema in Women Treated for Breast Cancer: A Pilot Study. Lymphatic research and biology, 18(1), 56-63. <https://doi.org/10.1089/lrb.2018.0067>

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8 cf. Dtsch Arztebl 2004; 101:A 2819-2824 [Heft 42]

9 cf. Anzarut, A., Olson, J., Singh, P., Rowe, B. H., & Tredget, E. E. (2009). The effectiveness of pressure garment therapy for the prevention of abnormal scarring after burn injury: a meta-analysis. *Journal of plastic, reconstructive & aesthetic surgery : JPRAS*, 62(1), 77-84. <https://doi.org/10.1016/j.bjps.2007.10.052>

10 cf. Li-Tsang, C. W., Zheng, Y. P., & Lau, J. C. (2010). A randomized clinical trial to study the effect of silicone gel dressing and pressure therapy on posttraumatic hypertrophic scars. *Journal of burn care & research : official publication of the American Burn Association*, 31(3), 448-457. <https://doi.org/10.1097/BCR.0b013e3181db52a7>

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©2022 Tous droits réservés Amoena Medizin-Orthopädie-Technik GmbH, Kapellenweg 36, 83064 Raubling, Allemagne
Amoena France S.A.S.